

兔抗 SF3A3 多克隆抗体

中文名称：兔抗 SF3A3 多克隆抗体

英文名称： Anti-SF3A3 rabbit polyclonal antibody

别名： splicing factor 3a subunit 3; PRP9; PRPF9; SAP61; SF3a60

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SF3A3

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes subunit 3 of the splicing factor 3a protein complex. The splicing factor 3a heterotrimer includes subunits 1, 2 and 3 and is necessary for the in vitro conversion of 15S U2 snRNP into an active 17S particle that performs pre-mRNA splicing. Subunit 3 interacts with subunit 1 through its amino-terminus while the zinc finger domain of subunit 3 plays a role in its binding to the 15S U2 snRNP. This gene has a pseudogene on chromosome 20. Alternative splicing results in multiple transcript variants.

Applications:	ELISA, WB, IHC
Name of antibody:	SF3A3
Immunogen:	Fusion protein of human SF3A3
Full name:	splicing factor 3a subunit 3
Synonyms:	PRP9; PRPF9; SAP61; SF3a60
SwissProt:	Q12874
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and Human liver cancer
IHC Recommend dilution:	30-150
WB Predicted band size:	59 kDa
WB Positive control:	Jurkat and Raji cell lysates
WB Recommended dilution:	500-2000



