

兔抗 PGD 多克隆抗体

中文名称：兔抗 PGD 多克隆抗体

英文名称： Anti-PGD rabbit polyclonal antibody

别名： 6PGD

抗原： PGD

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	6-phosphogluconate dehydrogenase is the second dehydrogenase in the pentose phosphate shunt. Deficiency of this enzyme is generally asymptomatic, and the inheritance of this disorder is autosomal dominant. Hemolysis results from combined deficiency of 6-phosphogluconate dehydrogenase and 6-phosphogluconolactonase suggesting a synergism of the two enzymopathies.
Applications:	WB
Name of antibody:	PGD
Immunogen:	Synthesized peptide derived from internal of human PGD.
Full name:	phosphogluconate dehydrogenase

Synonyms :	6PGD
SwissProt:	P52209
WB Predicted band size:	53 kDa
WB Positive control:	A549 cells lysate
WB Recommended dilution:	500-3000

