

兔抗 KAT7 多克隆抗体

中文名称:	兔抗 KAT7 多克隆抗体
英文名称:	Anti-KAT7 rabbit polyclonal antibody
别名:	lysine acetyltransferase 7; HBO1; HBOA; MYST2; ZC2HC7
相关类别:	一抗
储存:	冷冻 (-20℃)
宿主:	Rabbit
抗原:	KAT7
反应种属:	Human, Mouse, Rat
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is part of the multimeric HBO1 complex, which possesses histone H4-specific acetyltransferase activity. This activity is required for functional replication origins and is involved in transcriptional activation of some genes. In both cases, the acetylation of histone H4 helps unfold chromatin so that the DNA can be accessed and replicated or transcribed.
Applications:	ELISA, IHC
Name of antibody:	KAT7
Immunogen:	Synthetic peptide of human KAT7
Full name:	lysine acetyltransferase 7

Synonyms:	HBO1; HBOA; MYST2; ZC2HC7
SwissProt:	O95251
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	50-300

