

## 兔抗 SF3B6 多克隆抗体

中文名称: 兔抗 SF3B6 多克隆抗体

英文名称: Anti-SF3B6 rabbit polyclonal antibody

别名: P14; Ht006; SAP14; SAP14a; SF3B14; CGI-110; HSPC175; SF3B14a

储存: 冷冻 (-20℃) 避光

抗原: SF3B6

宿主: Rabbit

反应种属: Human Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a 14 kDa protein subunit of the splicing factor 3b complex. Splicing factor 3b associates with both the U2 and U11/U12 small nuclear ribonucleoprotein complexes (U2 snRNP) of spliceosomes. This 14 kDa protein interacts directly with subunit 1 of the splicing factor 3b complex. This 14 kDa protein also interacts directly with the adenosine that carries out the first transesterification step of splicing at the pre-mRNA branch site.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	SF3B6

<b>Immunogen:</b>	Synthesized peptide derived from C-terminal of human SF3B14.
<b>Full name:</b>	splicing factor 3b, subunit 6, 14kDa
<b>Synonyms :</b>	P14; Ht006; SAP14; SAP14a; SF3B14; CGI-110; HSPC175; SF3B14a
<b>SwissProt:</b>	Q9Y3B4
<b>WB Predicted band size:</b>	15 kDa
<b>WB Positive control:</b>	HepG2 cells lysate
<b>WB Recommended dilution:</b>	500-3000

