

## 兔抗 SFR1 多克隆抗体

中文名称：兔抗 SFR1 多克隆抗体

英文名称：Anti-SFR1 rabbit polyclonal antibody

别名：MEI5; MEIR5; C10orf78; bA373N18.1

相关类别：一抗

抗原：SFR1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Component of the SWI5-SFR1 complex, a complex required for double-strand break repair via homologous recombination (PubMed:21252223). Acts as a transcriptional modulator for ESR1 (PubMed:23874500).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SFR1
<b>Immunogen:</b>	Synthetic peptide of human SFR1
<b>Full name:</b>	SWI5 dependent homologous recombination repair protein 1
<b>Synonyms:</b>	MEI5; MEIR5; C10orf78; bA373N18.1
<b>SwissProt:</b>	Q86XK3

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	20-100

