

兔抗 WHRN 多克隆抗体

中文名称：兔抗 WHRN 多克隆抗体

英文名称： Anti-WHRN rabbit polyclonal antibody

别名： WI; CIP98; USH2D; DFNB31; PDZD7B

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： WHRN

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is thought to function in the organization and stabilization of stereocilia elongation and actin cytoskeletal assembly, based on studies of the related mouse gene. Mutations in this gene have been associated with autosomal recessive non-syndromic deafness and Usher Syndrome. Alternative splicing of this gene results in multiple transcript variants encoding different isoforms.
Applications:	ELISA, WB
Name of antibody:	WHRN
Immunogen:	Synthetic peptide of human WHRN

Full name:	whirlin
Synonyms:	WI; CIP98; USH2D; DFNB31; PDZD7B
SwissProt:	Q9P202
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	97 kDa
WB Positive control:	293T, PC-3 and TM4 cell lysates
WB Recommended dilution:	200-1000

