

## 兔抗 SFRP2 多克隆抗体

中文名称: 兔抗 SFRP2 多克隆抗体

英文名称: Anti-SFRP2 rabbit polyclonal antibody

别名: secreted frizzled related protein 2; FRP-2; SARP1; SDF-5

储存: 冷冻 (-20℃)

抗原: SFRP2

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. SFRPs act as soluble modulators of Wnt signaling. Methylation of this gene is a potential marker for the presence of colorectal cancer.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SFRP2
<b>Immunogen:</b>	Synthetic peptide of human SFRP2
<b>Full name:</b>	secreted frizzled related protein 2
<b>Synonyms:</b>	FRP-2; SARP1; SDF-5

<b>SwissProt:</b>	Q96HF1
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	25-100

