

兔抗 PGK2 多克隆抗体

中文名称：兔抗 PGK2 多克隆抗体

英文名称： Anti-PGK2 rabbit polyclonal antibody

别名： PGKB; PGKPS; dJ417L20.2

相关类别 一抗

抗原： PGK2

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is intronless, arose via retrotransposition of the phosphoglycerate kinase 1 gene, and is expressed specifically in the testis. Initially assumed to be a pseudogene, the encoded protein is actually a functional phosphoglycerate kinase that catalyzes the reversible conversion of 1,3-bisphosphoglycerate to 3-phosphoglycerate, during the Embden-Meyerhof-Parnas pathway of glycolysis, in the later stages of spermatogenesis.
Applications:	ELISA, WB, IHC
Name of antibody:	PGK2
Immunogen:	Fusion protein of human PGK2

Full name:	phosphoglycerate kinase 2
Synonyms :	PGKB; PGKPS; dJ417L20.2
SwissProt:	P07205
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human thyroid cancer and human liver cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	45 kDa
WB Positive control:	HepG2, 293T, hela and NIH/3T3 cells
WB Recommended dilution:	1000-5000



