

兔抗 SFRP5 多克隆抗体

中文名称：兔抗 SFRP5 多克隆抗体

英文名称：Anti-SFRP5 rabbit polyclonal antibody

别名：SARP3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SFRP5

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Secreted frizzled-related protein 5 (SFRP5) is a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. SFRPs act as soluble modulators of Wnt signaling. SFRP5 and SFRP1 may be involved in determining the polarity of photoreceptor cells in the retina. SFRP5 is highly expressed in the retinal pigment epithelium, and moderately expressed in the pancreas.
Applications:	ELISA, IHC
Name of antibody:	SFRP5

Immunogen:	Synthetic peptide of human SFRP5
Full name:	secreted frizzled-related protein 5
Synonyms:	SARP3
SwissProt:	Q5T4F7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and human thyroid cancer
IHC Recommend dilution:	25-100



