

## 兔抗 KATNAL1 多克隆抗体

中文名称：兔抗 KATNAL1 多克隆抗体

英文名称：Anti-KATNAL1 rabbit polyclonal antibody

别名：katanin catalytic subunit A1 like 1

相关类别：一抗 储存 冷冻 (-20℃)

宿主：Rabbit

抗原：KATNAL1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Regulates microtubule dynamics in Sertoli cells, a process that is essential for spermiogenesis and male fertility. Severs microtubules in an ATP-dependent manner, promoting rapid reorganization of cellular microtubule arrays (By similarity). Has microtubule-severing activity in vitro (PubMed:26929214).
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	KATNAL1
<b>Immunogen:</b>	Fusion protein of human KATNAL1
<b>Full name:</b>	katanin catalytic subunit A1 like 1
<b>SwissProt:</b>	Q9BW62
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human cervical cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	50-300
<b>WB Predicted band size:</b>	55 kDa
<b>WB Positive control:</b>	Mouse brain tissue lysate
<b>WB Recommended dilution:</b>	500-2000



