

## 兔抗 WIPF2 多克隆抗体

中文名称：兔抗 WIPF2 多克隆抗体

英文名称：Anti-WIPF2 rabbit polyclonal antibody

级别：WICH; WIRE

抗原：WIPF

储存：冷冻 (-20℃)

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a WASP interacting protein (WIP)-related protein. It has been shown that this protein has a role in the WASP-mediated organization of the actin cytoskeleton and that this protein is a potential link between the activated platelet-derived growth factor receptor and the actin polymerization machinery.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	WIPF2
<b>Immunogen:</b>	Synthetic peptide of human WIPF2
<b>Full name:</b>	WAS/WASL interacting protein family member 2

<b>Synonyms:</b>	WICH; WIRE
<b>SwissProt:</b>	Q8TF74
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and human lung cancer
<b>IHC Recommend dilution:</b>	10-50
<b>WB Predicted band size:</b>	46 kDa
<b>WB Positive control:</b>	Hela cell lysate
<b>WB Recommended dilution:</b>	500-2000



