

## 兔抗 SFT2D2 多克隆抗体

中文名称：兔抗 SFT2D2 多克隆抗体

英文名称：Anti-SFT2D2 rabbit polyclonal antibody

别名：UNQ512; dJ747L4.C1.2

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：SFT2D2

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex By similarity. UniProtKB P38166
<b>Applications:</b>	WB
<b>Name of antibody:</b>	SFT2D2
<b>Immunogen:</b>	Synthesized peptide derived from internal of human SFT2B.
<b>Full name:</b>	SFT2 domain containing 2
<b>Synonyms :</b>	UNQ512; dJ747L4.C1.2
<b>SwissProt:</b>	O95562

<b>WB Predicted band size:</b>	18 kDa
<b>WB Positive control:</b>	LOVO cells lysate
<b>WB Recommended dilution:</b>	500-3000

