

## 兔抗 WNK2 多克隆抗体

中文名称： 兔抗 WNK2 多克隆抗体

英文名称： Anti-WNK2 rabbit polyclonal antibody

别名： PRKWNK2; NY-CO-43; SDCCAG43; P/OKcl.13

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： WNK2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a cytoplasmic serine-threonine kinase that belongs to the protein kinase superfamily. The protein plays an important role in the regulation of electrolyte homeostasis, cell signaling survival, and proliferation. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	WNK2
<b>Immunogen:</b>	Synthetic peptide of human WNK2
<b>Full name:</b>	WNK lysine deficient protein kinase 2
<b>Synonyms :</b>	PRKWNK2; NY-CO-43; SDCCAG43; P/OKcl.13

<b>SwissProt:</b>	Q9Y3S1
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human thyroid cancer and Human prostate cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	243 kDa
<b>WB Positive control:</b>	A431 and MCF7 cells
<b>WB Recommended dilution:</b>	200-1000



