

兔抗 WNK4 多克隆抗体

中文名称：兔抗 WNK4 多克隆抗体

英文名称： Anti-WNK4 rabbit polyclonal antibody

别名： WNK lysine deficient protein kinase 4; PHA2B; PRKWINK4

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： WNK4

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the WNK family of serine-threonine protein kinases. The kinase is part of the tight junction complex in kidney cells, and regulates the balance between NaCl reabsorption and K(+) secretion. The kinase regulates the activities of several types of ion channels, cotransporters, and exchangers involved in electrolyte flux in epithelial cells. Mutations in this gene result in pseudohypoaldosteronism type IIB.
Applications:	ELISA, IHC
Name of antibody:	WNK4
Immunogen:	Synthetic peptide of human WNK4
Full name:	WNK lysine deficient protein kinase 4
Synonyms:	PHA2B; PRKWINK4

SwissProt:	Q96J92
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	25-100

