

## 兔抗 SGCA 多克隆抗体

中文名称：兔抗 SGCA 多克隆抗体

英文名称： Anti-SGCA rabbit polyclonal antibody

别名： sarcoglycan alpha; ADL; DAG2; 50DAG; DMDA2; LGMD2D; SCARMD1; adhalin

抗原： SGCA

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a component of the dystrophin-glycophorin complex (DGC), which is critical to the stability of muscle fiber membranes and to the linking of the actin cytoskeleton to the extracellular matrix. Its expression is thought to be restricted to striated muscle. Mutations in this gene result in type 2D autosomal recessive limb-girdle muscular dystrophy. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SGCA

<b>Immunogen:</b>	Fusion protein of human SGCA
<b>Full name:</b>	sarcoglycan alpha
<b>Synonyms:</b>	ADL; DAG2; 50DAG; DMDA2; LGMD2D; SCARMD1; adhalin
<b>SwissProt:</b>	Q16586
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and Human lung cancer
<b>IHC Recommend dilution:</b>	25-100



