

兔抗 WNT2B 多克隆抗体

中文名称： 兔抗 WNT2B 多克隆抗体

英文名称： Anti-WNT2B rabbit polyclonal antibody

别 名： WNT13

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： WNT2B

反应种属： Human, Mouse

标 记 物： Unconjugate

克隆类型： rabbit polyclonal

包装：

25 ul 400

100ul 1500

200 ul 2500

技术规格

Backg round:	This gene encodes a member of the wingless-type MMTV integration site (WNT) family of highly conserved, secreted signaling factors. WNT family members function in a variety of developmental processes including regul
-------------------------	---

	<p>ation of cell growth and differentiation and are characterized by a WNT-core domain. This gene may play a role in human development as well as carcinogenesis. Alternative splicing results in multiple transcript variants.</p>
Applications:	ELISA, IHC
Name of antibody:	WNT2B
Immunogen:	Synthetic peptide of human WNT2B
Full name:	wingless-type MMTV integration site family, member 2B
Synonyms:	WNT13
Swiss Prot:	Q93097
ELISA Recommended dilution:	1000-2000

IHC p ositive contro l:	Human esophagus cancer
IHC R ecom mend dilutio n:	10-50