

## 兔抗 PGRMC1 多克隆抗体

中文名称：兔抗 PGRMC1 多克隆抗体

英文名称：Anti-PGRMC1 rabbit polyclonal antibody

别名：MPR; HPR6.6

抗原：PGRMC1

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a putative membrane-associated progesterone steroid receptor. The protein is expressed predominantly in the liver and kidney.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PGRMC1
<b>Immunogen:</b>	Fusion protein of human PGRMC1
<b>Full name:</b>	progesterone receptor membrane component 1
<b>Synonyms:</b>	MPR; HPR6.6
<b>SwissProt:</b>	O00264
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human prostate cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	22 kDa
<b>WB Positive control:</b>	Human kidney tissue lysate
<b>WB Recommended dilution:</b>	200-1000



