

兔抗 WNT2 多克隆抗体

中文名称：兔抗 WNT2 多克隆抗体

英文名称：Anti-WNT2 rabbit polyclonal antibody

别名：IRP; INT1L1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：WNT2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the WNT gene family. The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. Alternatively spliced transcript variants have been identified for this gene.
Applications:	ELISA, IHC
Name of antibody:	WNT2
Immunogen:	Synthetic peptide of human WNT2
Full name:	wingless-type MMTV integration site family member 2

Synonyms :	IRP; INT1L1
SwissProt:	P09544
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human prostate cancer and Human thyroid cancer
IHC Recommend dilution:	25-100

