

兔抗 WNT5B 多克隆抗体

中文名称：兔抗 WNT5B 多克隆抗体

英文名称： Anti-WNT5B rabbit polyclonal antibody

别名： Wnt family member 5B

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： WNT5B

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It encodes a protein which shows 94% and 80% amino acid identity to the mouse Wnt5b protein and the human WNT5A protein, respectively. Alternative splicing of this gene generates 2 transcript variants.
Applications:	ELISA, IHC

Name of antibody:	WNT5B
Immunogen:	Synthetic peptide of human WNT5B
Full name:	Wnt family member 5B
SwissProt:	Q9H1J7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer
IHC Recommend dilution:	30-150

