

## 兔抗 PHB 多克隆抗体

中文名称：兔抗 PHB 多克隆抗体

英文名称： Anti-PHB rabbit polyclonal antibody

别名： PHB1

抗原： PHB

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is evolutionarily conserved, and its product is proposed to play a role in human cellular senescence and tumor suppression. Antiproliferative activity is reported to be localized to the 3' UTR, which is proposed to function as a trans-acting regulatory RNA. Several pseudogenes of this gene have been identified. Alternative splicing results in multiple transcript variants
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PHB
<b>Immunogen:</b>	Fusion protein of human PHB
<b>Full name:</b>	prohibitin

<b>Synonyms :</b>	PHB1
<b>SwissProt:</b>	P35232
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer and Human brain
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	30 kDa
<b>WB Positive control:</b>	293T and A549 cells, human hepatocellular carcinoma tissue and PC3 cells,hela cells and mouse heart tissue
<b>WB Recommended dilution:</b>	500-2000



