

兔抗 PHF1 多克隆抗体

中文名称：兔抗 PHF1 多克隆抗体

英文名称： Anti-PHF1 rabbit polyclonal antibody

别名： PCL1; PHF2; hPHF1; MTF2L2; TDRD19C

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： PHF1

宿主： Rabbit

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a Polycomb group protein. The protein is a component of a histone H3 lysine-27 (H3K27)-specific methyltransferase complex, and functions in transcriptional repression of homeotic genes. The protein is also recruited to double-strand breaks, and reduced protein levels results in X-ray sensitivity and increased homologous recombination. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	WB
Name of antibody:	PHF1

Immunogen:	Synthesized peptide derived from internal of human PHF1.
Full name:	PHD finger protein 1
Synonyms :	PCL1; PHF2; hPHF1; MTF2L2; TDRD19C
SwissProt:	O43189
WB Predicted band size:	62 kDa
WB Positive control:	Hela cells lysate
WB Recommended dilution:	500-3000

