

兔抗 WNT11 多克隆抗体

中文名称：兔抗 WNT11 多克隆抗体

英文名称：Anti-WNT11 rabbit polyclonal antibody

别名：HWNT11

相关类别：一抗

储存：冷冻（-20℃）

抗原：WNT11

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It encodes a protein which shows 97%, 85%, and 63% amino acid identity with mouse, chicken, and Xenopus Wnt11 protein, respectively. This gene may play roles in the development of skeleton, kidney and lung, and is considered to be a plausible candidate gene for High Bone Mass Syndrome

Applications:	ELISA, IHC
Name of antibody:	WNT11
Immunogen:	Fusion protein of human WNT11
Full name:	wingless-type MMTV integration site family, member 11
Synonyms :	HWNT11
SwissProt:	O96014
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human prostate cancer
IHC Recommend dilution:	25-100

