

兔抗 PHF3 多克隆抗体

中文名称: 兔抗 PHF3 多克隆抗体

英文名称: Anti-PHF3 rabbit polyclonal antibody

储 存: 冷冻 (-20℃) 避光

抗 原: PHF3

宿 主: Rabbit

反应种属: Human

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	This gene encodes a member of a PHD finger-containing gene family. This gene may function as a transcription factor and may be involved in glioblastomas development. Alternative splicing results in multiple transcript variants.
Applications:	WB, IHC
Name of antibody:	PHF3
Immunogen:	Synthesized peptide derived from C-terminal of human PHF3
Full name:	PHD finger protein 3
Synonyms :	0
SwissProt:	Q92576
IHC positive control:	Human breast carcinoma tissue

IHC Recommend dilution:	50-100
WB Predicted band size:	229 kDa
WB Positive control:	Hela cells lysate
WB Recommended dilution:	500-3000



