

兔抗 PHF5A 多克隆抗体

中文名称：兔抗 PHF5A 多克隆抗体

英文名称： Anti-PHF5A rabbit polyclonal antibody

别名： PHD finger protein 5A; INI; Rds3; SF3B7; SAP14b; SF3b14b; bK223H9.2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PHF5A

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a subunit of the splicing factor 3b protein complex. Splicing factor 3b, together with splicing factor 3a and a 12S RNA unit, forms the U2 small nuclear ribonucleoproteins complex (U2 snRNP). The splicing factor 3b/3a complex binds pre-mRNA upstream of the intron's branch site in a sequence-independent manner and may anchor the U2 snRNP to the pre-mRNA. The protein encoded by this gene contains a PHD-finger-like domain that is flanked by highly basic N- and C-termini. This protein belongs to the PHD-finger su

	perfamily and may act as a chromatin-associated protein. This gene has several pseudogenes on different chromosomes.
Applications:	ELISA, IHC
Name of antibody:	PHF5A
Immunogen:	Fusion protein of human PHF5A
Full name:	PHD finger protein 5A
Synonyms:	INI; Rds3; SF3B7; SAP14b; SF3b14b; bK223H9.2
SwissProt:	Q7RTV0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	100-200

