

兔抗 PHF6 多克隆抗体

中文名称：兔抗 PHF6 多克隆抗体

英文名称： Anti-PHF6 rabbit polyclonal antibody

别名： PHD finger protein 6; BFLS; BORJ; CENP-31

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PHF6

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene is a member of the plant homeodomain (PHD)-like finger (PHF) family. It encodes a protein with two PHD-type zinc finger domains, indicating a potential role in transcriptional regulation, that localizes to the nucleolus. Mutations affecting the coding region of this gene or the splicing of the transcript have been associated with Borjeson-Forssman-Lehmann syndrome (BFLS), a disorder characterized by cognitive disability, epilepsy, hypogonadism, hypometabolism, obesity, swelling of subcutaneous tissue of the face, narrow palpebral fissures,

	and large ears. Alternate splicing results in multiple transcript variants, encoding different isoforms.
Applications:	ELISA, IHC
Name of antibody:	PHF6
Immunogen:	Fusion protein of human PHF6
Full name:	PHD finger protein 6
Synonyms:	BFLS; BORJ; CENP-31
SwissProt:	Q8IWS0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer and Human lung cancer
IHC Recommend dilution:	100-300



