

兔抗 PHKB 多克隆抗体

中文名称：兔抗 PHKB 多克隆抗体

英文名称： Anti-PHKB rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： PHKB

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Phosphorylase kinase is a polymer of 16 subunits, four each of alpha, beta, gamma and delta. The alpha subunit includes the skeletal muscle and hepatic isoforms, encoded by two different genes. The beta subunit is the same in both the muscle and hepatic isoforms, encoded by this gene, which is a member of the phosphorylase b kinase regulatory subunit family. The gamma subunit also includes the skeletal muscle and hepatic isoforms, encoded by two different genes.
Applications:	ELISA, IHC
Name of antibody:	PHKB
Immunogen:	Fusion protein of human PHKB

Full name:	phosphorylase kinase, beta
SwissProt:	Q93100
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer and Human esophagus cancer
IHC Recommend dilution:	25-100

