

## 兔抗 PHKG1 多克隆抗体

中文名称: 兔抗 PHKG1 多克隆抗体

英文名称: Anti-PHKG1 rabbit polyclonal antibody

别名: PHKG

抗原: PHKG1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the Ser/Thr protein kinase family and encodes a protein with one protein kinase domain and two calmodulin-binding domains. This protein is the catalytic member of a 16 subunit protein kinase complex which contains equimolar ratios of 4 subunit types. The complex is a crucial glycogenolytic regulatory enzyme. This gene has two pseudogenes at chromosome 7q11.21 and one at chromosome 11p11.12.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PHKG1

<b>Immunogen:</b>	Synthetic peptide of human PHKG1
<b>Full name:</b>	phosphorylase kinase, gamma 1 (muscle)
<b>Synonyms:</b>	PHKG
<b>SwissProt:</b>	Q16816
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human brain and human liver cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	45 kDa
<b>WB Positive control:</b>	Human fetal muscle tissue
<b>WB Recommended dilution:</b>	200-1000



