

## 兔抗 PHLDA2 多克隆抗体

中文名称：兔抗 PHLDA2 多克隆抗体

英文名称： Anti-PHLDA2 rabbit polyclonal antibody

别名： IPL; BRW1C; BWR1C; HLDA2; TSSC3

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： PHLDA2

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene is located in a cluster of imprinted genes on chromosome 11p15.5, which is considered to be an important tumor suppressor gene region. Alterations in this region may be associated with the Beckwith-Wiedemann syndrome, Wilms tumor, rhabdomyosarcoma, adrenocortical carcinoma, and lung, ovarian, and breast cancer. This gene has been shown to be imprinted, with preferential expression from the maternal allele in placenta and liver.
<b>Applications:</b>	IHC
<b>Name of antibody:</b>	PHLDA2

<b>Immunogen:</b>	Synthesized peptide derived from internal of human PHLA2.
<b>Full name:</b>	pleckstrin homology-like domain, family A, member 2
<b>Synonyms :</b>	IPL; BRW1C; BWR1C; HLDA2; TSSC3
<b>SwissProt:</b>	Q53GA4
<b>IHC positive control:</b>	Human placenta tissue
<b>IHC Recommend dilution:</b>	50-100

