

兔抗 PHLDA3 多克隆抗体

中文名称：兔抗 PHLDA3 多克隆抗体

英文名称： Anti-PHLDA3 rabbit polyclonal antibody

别名： pleckstrin homology like domain family A member 3; TIH1

抗原： PHLDA3

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	p53/TP53-regulated repressor of Akt/AKT1 signaling. Represses AKT1 by preventing AKT1-binding to membrane lipids, thereby inhibiting AKT1 translocation to the cellular membrane and activation. Contributes to p53/TP53-dependent apoptosis by repressing AKT1 activity. Its direct transcription regulation by p53/TP53 may explain how p53/TP53 can negatively regulate AKT1. May act as a tumor suppressor.
Applications:	ELISA, IHC
Name of antibody:	PHLDA3

Immunogen:	Fusion protein of human PHLDA3
Full name:	pleckstrin homology like domain family A member 3
Synonyms:	TIH1
SwissProt:	Q9Y5J5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human thyroid cancer
IHC Recommend dilution:	100-300



