

兔抗 PHLPP1 多克隆抗体

中文名称：兔抗 PHLPP1 多克隆抗体

英文名称： Anti-PHLPP1 rabbit polyclonal antibody

别名： PHLPP1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PHLPP1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the serine/threonine phosphatase family. The encoded protein promotes apoptosis by dephosphorylating and inactivating the serine/threonine kinase Akt, and functions as a tumor suppressor in multiple types of cancer. Increased expression of this gene may also play a role in obesity and type 2 diabetes by interfering with Akt-mediated insulin signaling.
Applications:	ELISA, IHC
Name of antibody:	PHLPP1
Immunogen:	Fusion protein of human PHLPP1

Full name:	PH domain and leucine rich repeat protein phosphatase 1
Synonyms :	PHLPP1
SwissProt:	O60346
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human esophagus cancer and human brain
IHC Recommend dilution:	25-100

