

兔抗 PHPT1 多克隆抗体

中文名称：兔抗 PHPT1 多克隆抗体

英文名称： Anti-PHPT1 rabbit polyclonal antibody

别名： phosphohistidine phosphatase 1; PHP; PHP14; CGI-202; HSPC141; HEL-S-132P

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PHPT1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an enzyme that catalyzes the reversible dephosphorylation of histidine residues in proteins. It may be involved in the dephosphorylation of G-beta and ATP citrate lyase and in negatively regulating CD4 T lymphocytes by dephosphorylation and inhibition of KCa3.1 channels. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	PHPT1
Immunogen:	Fusion protein of human PHPT1
Full name:	phosphohistidine phosphatase 1

Synonyms:	PHP; PHP14; CGI-202; HSPC141; HEL-S-132P
SwissProt:	Q9NRX4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	200-400
WB Predicted band size:	14 kDa
WB Positive control:	293T cell lysate
WB Recommended dilution:	1000-5000



