

兔抗 PIGQ 多克隆抗体

中文名称：兔抗 PIGQ 多克隆抗体

英文名称：Anti-PIGQ rabbit polyclonal antibody

别名：phosphatidylinositol glycan anchor biosynthesis class Q; GPI1; c407A10.1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PIGQ

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is involved in the first step in glycosylphosphatidylinositol (GPI)-anchor biosynthesis. The GPI-anchor is a glycolipid found on many blood cells and serves to anchor proteins to the cell surface. This gene encodes a N-acetylglucosaminyl transferase component that is part of the complex that catalyzes transfer of N-acetylglucosamine (GlcNAc) from UDP-GlcNAc to phosphatidylinositol (PI). Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC

Name of antibody:	PIGQ
Immunogen:	Fusion protein of human PIGQ
Full name:	phosphatidylinositol glycan anchor biosynthesis class Q
Synonyms:	GPI1; c407A10.1
SwissProt:	Q9BRB3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	50-200

