

## 兔抗 PIK3CD 多克隆抗体

中文名称：兔抗 PIK3CD 多克隆抗体

英文名称： Anti-PIK3CD rabbit polyclonal antibody

别名： APDS; PI3K; IMD14; p110D; P110DELTA

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Rabbit

抗原： PIK3CD

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Phosphoinositide 3-kinases (PI3Ks) phosphorylate inositol lipids and are involved in the immune response. The protein encoded by this gene is a class I PI3K found primarily in leukocytes. Like other class I PI3Ks (p110-alpha p110-beta, and p110-gamma), the encoded protein binds p85 adapter proteins and GTP-bound RAS. However, unlike the other class I PI3Ks, this protein phosphorylates itself, not p85 protein.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PIK3CD
<b>Immunogen:</b>	Synthetic peptide of human PIK3CD

<b>Full name:</b>	phosphatidylinositol-4,5-bisphosphate 3-kinase, catalytic subunit delta
<b>Synonyms :</b>	APDS; PI3K; IMD14; p110D; P110DELTA
<b>SwissProt:</b>	O00329
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human cervical cancer
<b>IHC Recommend dilution:</b>	25-100

