

兔抗 PIK3R4 多克隆抗体

中文名称：兔抗 PIK3R4 多克隆抗体

英文名称：Anti-PIK3R4 rabbit polyclonal antibody

别名：p150; VPS15

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PIK3R4

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Phosphoinositide 3-kinase regulatory subunit 4, also known as PI3-kinase regulatory subunit 4 or PI3-kinase p150 subunit or phosphoinositide 3-kinase adaptor protein, or VPS15 is an enzyme that in humans is encoded by the PIK3R4 gene. The PI3-kinases regulate cellular signaling networks that are involved in processes linked to the survival, growth, proliferation, metabolism and specialized differentiated functions of cells. The subversion of this network is common in cancer and has also been linked to disorders of inflammation.
Applications:	ELISA, IHC
Name of antibody:	PIK3R4
Immunogen:	Synthetic peptide of human PIK3R4

Full name:	phosphoinositide-3-kinase, regulatory subunit 4
Synonyms :	p150; VPS15
SwissProt:	Q99570
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human lymphoma and Human thyroid cancer
IHC Recommend dilution:	100-300

