

## 兔抗 WSB2 多克隆抗体

中文名称：兔抗 WSB2 多克隆抗体

英文名称：Anti-WSB2 rabbit polyclonal antibody

别名：SBA2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：WSB2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the WD-protein subfamily. The encoded protein contains five WD-repeats spanning most of the protein and an SOCS box in the C-terminus. The SOCS box may act as a bridge between specific substrate-binding domains and E3 ubiquitin protein ligases. Three transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	WSB2
<b>Immunogen:</b>	Fusion protein of human WSB2
<b>Full name:</b>	WD repeat and SOCS box containing 2

<b>Synonyms:</b>	SBA2
<b>SwissProt:</b>	Q9NYS7
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	45 kDa
<b>WB Positive control:</b>	PC-3 cell lysate
<b>WB Recommended dilution:</b>	200-1000

