

兔抗 SGPL1 多克隆抗体

中文名称：兔抗 SGPL1 多克隆抗体

英文名称： Anti-SGPL1 rabbit polyclonal antibody

别名： SPL; S1PL

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SGPL1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Sphingosine-1-phosphate lyase 1 is an enzyme that in humans is encoded by the SGPL1 gene. Cleaves phosphorylated sphingoid bases (PSBs), such as sphingosine-1-phosphate, into fatty aldehydes and phosphoethanolamine. Elevates stress-induced ceramide production and apoptosis.
Applications:	ELISA, WB, IHC
Name of antibody:	SGPL1
Immunogen:	Synthetic peptide of human SGPL1
Full name:	sphingosine-1-phosphate lyase 1
Synonyms :	SPL; S1PL
SwissProt:	O95470

ELISA Recommended dilution:	1000-2000
IHC positive control:	Human breast cancer and Human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	64 kDa
WB Positive control:	Human placenta tissue
WB Recommended dilution:	200-1000



