

## 兔抗 SGTA 多克隆抗体

中文名称：兔抗 SGTA 多克隆抗体

英文名称： Anti-SGTA rabbit polyclonal antibody

别名： small glutamine rich tetratricopeptide repeat containing alpha; SGT; hSGT; alphaSGT

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SGTA

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein which is capable of interacting with the major nonstructural protein of parvovirus H-1 and 70-kDa heat shock cognate protein; however, its function is not known. Since this transcript is expressed ubiquitously in various tissues, this protein may serve a housekeeping function.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	SGTA
<b>Immunogen:</b>	Fusion protein of human SGTA
<b>Full name:</b>	small glutamine rich tetratricopeptide repeat containin

	g alpha
<b>Synonyms:</b>	SGT; hSGT; alphaSGT
<b>SwissProt:</b>	O43765
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human gastric cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	34 kDa
<b>WB Positive control:</b>	LNCAP and A172 cell lysates
<b>WB Recommended dilution:</b>	1000-5000



