

兔抗 WTAP 多克隆抗体

中文名称：兔抗 WTAP 多克隆抗体

英文名称：Anti-WTAP rabbit polyclonal antibody

别名：Mum2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：WTAP

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The Wilms tumor suppressor gene WT1 appears to play a role in both transcriptional and posttranscriptional regulation of certain cellular genes. This gene encodes a WT1-associated protein, which is a ubiquitously expressed nuclear protein. Like WT1 protein, this protein is localized throughout the nucleoplasm as well as in speckles and partially colocalizes with splicing factors. Alternative splicing of this gene results in several transcript variants encoding three different isoforms.
Applications:	ELISA, IHC
Name of antibody:	WTAP

Immunogen:	Full length fusion protein
Full name:	Wilms tumor 1 associated protein
Synonyms:	Mum2
SwissProt:	Q15007
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human brain
IHC Recommend dilution:	25-100



