

兔抗 SH2B2 多克隆抗体

中文名称：兔抗 SH2B2 多克隆抗体

英文名称：Anti-SH2B2 rabbit polyclonal antibody

别名：APS

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SH2B2

反应种属：Human, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is expressed in B lymphocytes and contains pleckstrin homology and src homology 2 (SH2) domains. In Burkitt's lymphoma cell lines, it is tyrosine-phosphorylated in response to B cell receptor stimulation. Because it binds Shc independent of stimulation and Grb2 after stimulation, it appears to play a role in signal transduction from the receptor to the Shc/Grb2 pathway.
Applications:	ELISA, IHC
Name of antibody:	SH2B2
Immunogen:	Synthetic peptide of human SH2B2

Full name:	SH2B adaptor protein 2
Synonyms :	APS
SwissProt:	O14492
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and Human gastric cancer
IHC Recommend dilution:	15-50

