

兔抗 SH2D1B 多克隆抗体

中文名称：兔抗 SH2D1B 多克隆抗体

英文名称： Anti-SH2D1B rabbit polyclonal antibody

别名： EAT2

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SH2D1B

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Plays a role in controlling signal transduction through at least four receptors, CD84, SLAMF1, LY9 and CD244, expressed on the surface of professional antigen-presenting cells.
Applications:	IHC
Name of antibody:	SH2D1B
Immunogen:	Synthesized peptide derived from internal of human SH2D1B.
Full name:	SH2 domain containing 1B
Synonyms :	EAT2
SwissProt:	O14796

IHC positive control:	Human lung carcinoma tissue
IHC Recommend dilution:	50-100

