

兔抗 WWP1 多克隆抗体

中文名称: 兔抗 WWP1 多克隆抗体

英文名称: Anti-WWP1 rabbit polyclonal antibody

别名: AIP5; Tiul1; hSDRP1

抗原: WWP1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

WW domain-containing proteins are found in all eukaryotes and play an important role in the regulation of a wide variety of cellular functions such as protein degradation, transcription, and RNA splicing. This gene encodes a protein which contains 4 tandem WW domains and a HECT (homologous to the E6-associated protein carboxyl terminus) domain. The encoded protein belongs to a family of NEDD4-like proteins, which are E3 ubiquitin-ligase molecules and regulate key trafficking decisions, including targeting of proteins to proteosomes or lysosomes.

Applications:	ELISA, WB, IHC
Name of antibody:	WWP1
Immunogen:	Fusion protein of human WWP1
Full name:	WW domain containing E3 ubiquitin protein ligase 1
Synonyms :	AIP5; Tiu1; hSDRP1
SwissProt:	Q9H0M0
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	30-150
WB Predicted band size:	105 kDa
WB Positive control:	231 and hela cell
WB Recommended dilution:	500-2000



