

兔抗 SH3GLB1 多克隆抗体

中文名称：兔抗 SH3GLB1 多克隆抗体

英文名称： Anti-SH3GLB1 rabbit polyclonal antibody

别名： Bif-1; CGI-61; PPP1R70; dJ612B15.2

相关类别： 一抗

抗原： SH3GLB1

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a SRC homology 3 domain-containing protein. The encoded protein interacts with the proapoptotic member of the Bcl-2 family, Bcl-2-associated X protein (Bax) and may be involved in regulating apoptotic signaling pathways. This protein may also be involved in maintaining mitochondrial morphology. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	SH3GLB1
Immunogen:	Fusion protein of human SH3GLB1
Full name:	SH3-domain GRB2-like endophilin B1

Synonyms:	Bif-1; CGI-61; PPP1R70; dJ612B15.2
SwissProt:	Q9Y371
IHC positive control:	Human gastric cancer and human thyroid cancer
IHC Recommend dilution:	15-50
WB Predicted band size:	41 kDa
WB Positive control:	Mouse heart tissue lysate
WB Recommended dilution:	200-1000



