

兔抗 XAF1 多克隆抗体

中文名称: 兔抗 XAF1 多克隆抗体

英文名称: Anti-XAF1 rabbit polyclonal antibody

别名: BIRC4BP; XIAPAF1; HSXIAPAF1

储存: 冷冻 (-20℃)

抗原: XAF1

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a protein which binds to and counteracts the inhibitory effect of a member of the IAP (inhibitor of apoptosis) protein family. IAP proteins bind to and inhibit caspases which are activated during apoptosis. The proportion of IAPs and proteins which interfere with their activity, such as the encoded protein, affect the progress of the apoptosis signaling pathway. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	XAF1

Immunogen:	Synthetic peptide of human XAF1
Full name:	XIAP associated factor 1
Synonyms:	BIRC4BP; XIAPAF1; HSXIAPAF1
SwissProt:	Q6GPH4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human prostate cancer and human tonsil
IHC Recommend dilution:	10-50



