

兔抗 XBP1 多克隆抗体

中文名称：兔抗 XBP1 多克隆抗体

英文名称： Anti-XBP1 rabbit polyclonal antibody

别名： XBP2; TREB5; XBP-1

抗原： XBP1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a transcription factor that regulates MHC class II genes by binding to a promoter element referred to as an X box. This gene product is a bZIP protein, which was also identified as a cellular transcription factor that binds to an enhancer in the promoter of the T cell leukemia virus type 1 promoter. It may increase expression of viral proteins by acting as the DNA binding partner of a viral transactivator. It has been found that upon accumulation of unfolded proteins in the endoplasmic reticulum (ER), the mRNA of this gene is processed to an active form by an unconventional splicing mechanism that is mediated by the endonuclease inositol-requiring enzyme

	1 (IRE1). The resulting loss of 26 nt from the spliced mRNA causes a frame-shift and an isoform XBP1(S), which is the functionally active transcription factor.
Applications:	ELISA, IHC
Name of antibody:	XBP1
Immunogen:	Fusion protein of human XBP1
Full name:	X-box binding protein 1
Synonyms :	XBP2; TREB5; XBP-1
SwissProt:	P17861
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human breast cancer and human thyroid cancer
IHC Recommend dilution:	50-200

