

## 兔抗 PIP5K1C 多克隆抗体

中文名称：兔抗 PIP5K1C 多克隆抗体

英文名称：Anti-PIP5K1C rabbit polyclonal antibody

别名：LCCS3; PIPKIg\_v4; PIP5Kgamma; PIP5K-GAMMA

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PIP5K1C

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This locus encodes a type I phosphatidylinositol 4-phosphate 5-kinase. The encoded protein catalyzes phosphorylation of phosphatidylinositol 4-phosphate, producing phosphatidylinositol 4,5-bisphosphate. This enzyme is found at synapses and has been found to play roles in endocytosis and cell migration. Mutations at this locus have been associated with lethal congenital contractural syndrome. Alternatively spliced transcript variants encoding different isoforms have been described.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PIP5K1C

<b>Immunogen:</b>	Synthetic peptide of human PIP5K1C
<b>Full name:</b>	phosphatidylinositol-4-phosphate 5-kinase, type I, gamma
<b>Synonyms :</b>	LCCS3; PIPKIg_v4; PIP5Kgamma; PIP5K-GAMMA
<b>SwissProt:</b>	O60331
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer and Human brain
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	73 kDa
<b>WB Positive control:</b>	A431 and hela cells
<b>WB Recommended dilution:</b>	500-2000



