

兔抗 PIPOX 多克隆抗体

中文名称: 兔抗 PIPOX 多克隆抗体

英文名称: Anti-PIPOX rabbit polyclonal antibody

别名: LPIPOX

储存: 冷冻 (-20℃) 避光

抗原: PIPOX

宿主: Rabbit

反应种属: Human Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

| | |
|---------------------------------|---|
| Background: | Metabolizes sarcosine, L-pipecolic acid and L-proline. |
| Applications: | WB |
| Name of antibody: | PIPOX |
| Immunogen: | Synthesized peptide derived from internal of human PIPOX. |
| Full name: | pipecolic acid oxidase |
| Synonyms : | LPIPOX |
| SwissProt: | Q9P0Z9 |
| WB Predicted band size: | 44 kDa |
| WB Positive control: | NIH/3T3 cells and RAW264.7 cells lysates |
| WB Recommended dilution: | 500-3000 |

