

兔抗 PIP 多克隆抗体

中文名称：兔抗 PIP 多克隆抗体

英文名称： Anti-PIP rabbit polyclonal antibody

别名： GPIP4; GCDFP15; GCDFP-15

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： PIP

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

| | |
|------------------------------------|---|
| Background: | The hormonally responsive prolactin-inducible protein (PIP) gene is expressed in benign and malignant breast tumor tissues and in some normal exocrine organs such as sweat, salivary, and lacrimal glands. |
| Applications: | ELISA, IHC |
| Name of antibody: | PIP |
| Immunogen: | Fusion protein of human PIP |
| Full name: | prolactin-induced protein |
| Synonyms : | GPIP4; GCDFP15; GCDFP-15 |
| SwissProt: | P12273 |
| ELISA Recommended dilution: | 2000-10000 |

| | |
|--------------------------------|---|
| IHC positive control: | Human breast cancer and Human cervical cancer |
| IHC Recommend dilution: | 100-300 |

