

## 兔抗 SHC3 多克隆抗体

中文名称：兔抗 SHC3 多克隆抗体

英文名称： Anti-SHC3 rabbit polyclonal antibody

别名： SHC adaptor protein 3; RAI; NSHC; SHCC; N-Shc

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SHC3

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Signaling adapter that couples activated growth factor receptors to signaling pathway in neurons. Involved in the signal transduction pathways of neurotrophin-activated Trk receptors in cortical neurons.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	SHC3
<b>Immunogen:</b>	Fusion protein of human SHC3
<b>Full name:</b>	SHC adaptor protein 3
<b>Synonyms:</b>	RAI; NSHC; SHCC; N-Shc
<b>SwissProt:</b>	Q92529
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	64 kDa
<b>WB Positive control:</b>	PC3 cell lysate
<b>WB Recommended dilution:</b>	500-2000



