

兔抗 PITX1 多克隆抗体

中文名称：兔抗 PITX1 多克隆抗体

英文名称： Anti-PITX1 rabbit polyclonal antibody

别名： BFT; CCF; POTX; PTX1; LBNBG

相关类别： 一抗

抗原： PITX1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a member of the RIEG/PITX homeobox family, which is in the bicoid class of homeodomain proteins. Members of this family are involved in organ development and left-right asymmetry. This protein acts as a transcriptional regulator involved in basal and hormone-regulated activity of prolactin.
Applications:	WB
Name of antibody:	PITX1
Immunogen:	Synthesized peptide derived from internal of human PITX1.
Full name:	paired-like homeodomain 1
Synonyms :	BFT; CCF; POTX; PTX1; LBNBG
SwissProt:	P78337

WB Predicted band size:	34 kDa
WB Positive control:	Jurkat cells and COLO cells lysates
WB Recommended dilution:	500-3000

