

## 兔抗 PITX3 多克隆抗体

中文名称：兔抗 PITX3 多克隆抗体

英文名称：Anti-PITX3 rabbit polyclonal antibody

别名：PTX3; CTPP4; CTRCT11

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PITX3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the RIEG/PITX homeobox family, which is in the bicoid class of homeodomain proteins. Members of this family act as transcription factors. This protein is involved in lens formation during eye development. Mutations of this gene have been associated with anterior segment mesenchymal dysgenesis and congenital cataracts.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PITX3
<b>Immunogen:</b>	Synthetic peptide of human PITX3
<b>Full name:</b>	paired-like homeodomain 3
<b>Synonyms :</b>	PTX3; CTPP4; CTRCT11

<b>SwissProt:</b>	O75364
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human lymphoma and Human brain
<b>IHC Recommend dilution:</b>	50-200

